

ATF5

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:31kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human ATF5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ATFX; HMFN0395;

To place an order, please [Click HERE](#).

Catalog #:A6859

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:22809

Isotype:IgG

Swiss Prot:Q9Y2D1

Purity:Affinity purification

For research use only.