

SREBF2

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:200 - 1:1000 IHC 1:50 - 1:100

Calculated MW:124kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human SREBF2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

SREBP2; bHLHd2; SREBP-2;

Catalog #:A6860

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6721

Isotype:IgG

Swiss Prot:Q12772

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the a ubiquitously expressed transcription factor that controls cholesterol homeostasis by regulating transcription of sterol-regulated genes. The encoded protein contains a basic helix-loop-helix-leucine zipper (bHLH-Zip) domain and binds the sterol regulatory element 1 motif. Alternate splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).