

BSG

Reactivity:Human

Tested applications:WB IHC ICC IF

Recommended Dilution:WB 1:500 - 1:1000 IHC 1:50 - 1:200 ICC 1:100 - 1:1000 IF 1:50 - 1:100

Calculated MW:42kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human BSG

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

M6; OK; 5F7; TCSF; CD147; EMMPRIN

Catalog #:A6861

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:682

Isotype:IgG

Swiss Prot:P35613

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a plasma membrane protein that is important in spermatogenesis, embryo implantation, neural network formation, and tumor progression. The encoded protein is also a member of the immunoglobulin superfamily. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).