

BMP3

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW:53kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human BMP3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BMP-3A;

Background:

BMP3 belongs to the transforming growth factor-beta (TGFB) superfamily. Bone morphogenic protein, also known as osteogenin, induces bone formation.

To place an order, please [Click HERE](#).

Catalog #:A6877

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:651

Isotype:IgG

Swiss Prot:P12645

Purity:Affinity purification

For research use only.