

## CKS1B

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**10kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human CKS1B

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

CKS1; ckshs1; PNAS-16; PNAS-18;

**Background:**

CKS1B protein binds to the catalytic subunit of the cyclin dependent kinases and is essential for their biological function. The CKS1B mRNA is found to be expressed in different patterns through the cell cycle in HeLa cells, which reflects a specialized role for the encoded protein. At least two transcript variants have been identified for this gene, and it appears that only one of them encodes a protein.

*To place an order, please [Click HERE](#).*

**Catalog #:**A6885

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**1163

**Isotype:**IgG

**Swiss Prot:**P61024

**Purity:**Affinity purification

For research use only.