

## DACH1

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**79kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human DACH1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

DACH;

**Background:**

This gene encodes a chromatin-associated protein that associates with other DNA-binding transcription factors to regulate gene expression and cell fate determination during development. The protein contains a Ski domain that is highly conserved from Drosophila to human. Expression of this gene is lost in some forms of metastatic cancer, and is correlated with poor prognosis. Multiple transcript variants encoding different isoforms have been found for this gene.

*To place an order, please [Click HERE](#).*

**Catalog #:**A6895

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**1602

**Isotype:**IgG

**Swiss Prot:**Q9UI36

**Purity:**Affinity purification

For research use only.