

DHODH

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW:43kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human DHODH

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

URA1; POADS; DHODEhase;

Catalog #:A6899

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1723

Isotype:IgG

Swiss Prot:Q02127

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene catalyzes the fourth enzymatic step, the ubiquinone-mediated oxidation of dihydroorotate to orotate, in de novo pyrimidine biosynthesis. This protein is a mitochondrial protein located on the outer surface of the inner mitochondrial membrane.

To place an order, please [Click HERE](#).