

DLST

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:49kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human DLST

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DLTS;

Catalog #:A6901

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1743

Isotype:IgG

Swiss Prot:P36957

Purity:Affinity purification

For research use only.

Background:

This gene encodes a mitochondrial protein that belongs to the 2-oxoacid dehydrogenase family.

This protein is one of the three components (the E2 component) of the 2-oxoglutarate dehydrogenase complex that catalyzes the overall conversion of 2-oxoglutarate to succinyl-CoA and CO(2). Alternatively spliced transcript variants have been found for this gene.

To place an order, please [Click HERE](#).