

## ELF1

---

**Reactivity:**Human Mouse

**Tested applications:**WB IHC IF

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

**Calculated MW:**67kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human ELF1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Catalog #:**A6904

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**1997

**Isotype:**IgG

**Swiss Prot:**P32519

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes an E26 transformation-specific related transcription factor. The encoded protein is primarily expressed in lymphoid cells and acts as both an enhancer and a repressor to regulate transcription of various genes. Alternative splicing results in multiple transcript variants.

**To place an order, please [Click HERE](#).**