

ENPEP

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:109kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human ENPEP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

APA; CD249; gp160;

To place an order, please [Click HERE](#).

Catalog #:A6905

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2028

Isotype:IgG

Swiss Prot:Q07075

Purity:Affinity purification

For research use only.