

## GIF

---

**Reactivity:**Human Mouse

**Tested applications:**WB IHC IF

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

**Calculated MW:**45kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human GIF

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

IF; INF; IFMH; TCN3;

**Catalog #:**A6914

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**2694

**Isotype:**IgG

**Swiss Prot:**P27352

**Purity:**Affinity purification

For research use only.

**Background:**

This gene is a member of the cobalamin transport protein family. It encodes a glycoprotein secreted by parietal cells of the gastric mucosa and is required for adequate absorption of vitamin B12. Vitamin B12 is necessary for erythrocyte maturation and mutations in this gene may lead to congenital pernicious anemia.

*To place an order, please [Click HERE](#).*