

GPS1

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:56kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human GPS1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CSN1; COPS1;

Catalog #:A6917

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2873

Isotype:IgG

Swiss Prot:Q13098

Purity:Affinity purification

For research use only.

Background:

This gene is known to suppress G-protein and mitogen-activated signal transduction in mammalian cells. The encoded protein shares significant similarity with Arabidopsis FUS6, which is a regulator of light-mediated signal transduction in plant cells. Two alternatively spliced transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).