

HABP2

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:62kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human HABP2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FSAP; HABP; PHBP; HGFAL;

Catalog #:A6919

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:3026

Isotype:IgG

Swiss Prot:Q14520

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is an extracellular serine protease that binds hyaluronic acid and is involved in cell adhesion. The encoded protein is synthesized as a single chain, but then undergoes an autoproteolytic event to form the functional heterodimer. Further autoproteolysis leads to smaller, inactive peptides. This protease is known to cleave urinary plasminogen activator, coagulation factor VII, and the alpha and beta chains of fibrinogen, but not prothrombin, plasminogen, or the gamma chain of fibrinogen. Two transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).