

MYBPC2

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:128kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human MYBPC2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MYBPC; MYBPCF;

Background:

This gene encodes a member of the myosin-binding protein C family. This family includes the fast-, slow- and cardiac-type isoforms, each of which is a myosin-associated protein found in the cross-bridge-bearing zone (C region) of A bands in striated muscle. The protein encoded by this locus is referred to as the fast-type isoform. Mutations in the related but distinct genes encoding the slow-type and cardiac-type isoforms have been associated with distal arthrogyrosis, type 1 and hypertrophic cardiomyopathy, respectively.

To place an order, please [Click HERE](#).

Catalog #:A6934

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4606

Isotype:IgG

Swiss Prot:Q14324

Purity:Affinity purification

For research use only.