

RBBP6

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 202kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RBBP6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PACT; MY038; P2P-R; RBQ-1; SNAMA;

Catalog #: A6966

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5930

Isotype: IgG

Swiss Prot: Q7Z6E9

Purity: Affinity purification

For research use only.

Background:

The retinoblastoma tumor suppressor (pRB) protein binds with many other proteins. In various human cancers, pRB suppresses cellular proliferation and is inactivated. Cell cycle-dependent phosphorylation regulates the activity of pRB. This gene encodes a protein which binds to underphosphorylated but not phosphorylated pRB. Multiple alternatively spliced transcript variants that encode different isoforms have been found for this gene.

To place an order, please [Click HERE](#).