

SGCB

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW:35kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human SGCB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

A3b; SGC; LGMD2E;

Background:

This gene encodes a member of the sarcoglycan family. Sarcoglycans are transmembrane components in the dystrophin-glycoprotein complex which help stabilize the muscle fiber membranes and link the muscle cytoskeleton to the extracellular matrix. Mutations in this gene have been associated with limb-girdle muscular dystrophy.

To place an order, please [Click HERE](#).

Catalog #:A6979

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6443

Isotype:IgG

Swiss Prot:Q16585

Purity:Affinity purification

For research use only.