

SNRPD2

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:14kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human SNRPD2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

SMD2; Sm-D2; SNRPD1;

Catalog #:A6983

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6633

Isotype:IgG

Swiss Prot:P62316

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene belongs to the small nuclear ribonucleoprotein core protein family. It is required for pre-mRNA splicing and small nuclear ribonucleoprotein biogenesis. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).