

SSBP1

Reactivity:Human

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:17kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human SSBP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SSBP; mtSSB; Mt-SSB; SOSS-B1;

Background:

SSBP1 is a housekeeping gene involved in mitochondrial biogenesis (Tiranti et al., 1995)

To place an order, please [Click HERE](#).

Catalog #:A6987

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6742

Isotype:IgG

Swiss Prot:Q04837

Purity:Affinity purification

For research use only.