

## STAT4

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC IF

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

**Calculated MW:**86kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human STAT4

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

ko

**Synonym:**

SLEB11;

**Catalog #:**A6991

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**6775

**Isotype:**IgG

**Swiss Prot:**Q14765

**Purity:**Affinity purification

For research use only.

**Background:**

The protein encoded by this gene is a member of the STAT family of transcription factors. In response to cytokines and growth factors, STAT family members are phosphorylated by the receptor associated kinases, and then form homo- or heterodimers that translocate to the cell nucleus where they act as transcription activators. This protein is essential for mediating responses to IL12 in lymphocytes, and regulating the differentiation of T helper cells. Mutations in this gene may be associated with systemic lupus erythematosus and rheumatoid arthritis. Alternate splicing results in multiple transcript variants that encode the same protein.

**To place an order, please [Click HERE](#).**