

COPS3

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:48kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human COPS3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CSN3; SGN3;

Catalog #:A7017

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8533

Isotype:IgG

Swiss Prot:Q9UNS2

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene possesses kinase activity that phosphorylates regulators involved in signal transduction. It phosphorylates I kappa-Balpha, p105, and c-Jun. It acts as a docking site for complex-mediated phosphorylation. The gene is located within the Smith-Magenis syndrome region on chromosome 17. Two transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).