

## NPFF

---

**Reactivity:**Human

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**12kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human NPFF

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

FMRFAL;

**Background:**

FMFRamide-related protein precursor plays a role in the regulation of heart rate and blood pressure and the modulation of morphine-induced antinociception. FMRFAL encodes a preproprotein which is cleaved to form two active peptides with similar function.

*To place an order, please [Click HERE](#).*

**Catalog #:**A7020

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**8620

**Isotype:**IgG

**Swiss Prot:**O15130

**Purity:**Affinity purification

For research use only.