

## EIF2B2

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**39kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human EIF2B2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

EIF2B; EIF-2Bbeta;

**Background:**

This gene encodes the beta subunit of eukaryotic initiation factor-2B (EIF2B). EIF2B is involved in protein synthesis and exchanges GDP and GTP for its activation and deactivation.

*To place an order, please [Click HERE](#).*

**Catalog #:**A7027

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**8892

**Isotype:**IgG

**Swiss Prot:**P49770

**Purity:**Affinity purification

For research use only.