

TNFSF18

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50- 1:200

Calculated MW:23kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human TNFSF18

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AITRL;GITRL;TL6; Tumor necrosis factor ligand superfamily member 18

Background:

Cytokine that binds to TNFRSF18/AITR/GITR. Important for interactions between activated T-lymphocytes and endothelial cells and may modulate T-lymphocyte survival in peripheral tissues.

To place an order, please [Click HERE](#).

Catalog #:A7028

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8995

Isotype:IgG

Swiss Prot:Q9UNG2

Purity:Affinity purification

For research use only.