

USP8

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 128kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human USP8

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

UBPY; SPG59; HumORF8;

Catalog #: A7031

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9101

Isotype: IgG

Swiss Prot: P40818

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein that belongs to the ubiquitin-specific processing protease family of proteins. The encoded protein is thought to regulate the morphology of the endosome by ubiquitination of proteins on this organelle and is involved in cargo sorting and membrane trafficking at the early endosome stage. This protein is required for the cell to enter the S phase of the cell cycle and also functions as a positive regulator in the Hedgehog signaling pathway in development. Pseudogenes of this gene are present on chromosomes 2 and 6. Alternate splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).