

COPB2

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:102kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human COPB2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

beta¹-COP;

Background:

The Golgi coatomer complex (see MIM 601924) constitutes the coat of nonclathrin-coated vesicles and is essential for Golgi budding and vesicular trafficking. It consists of 7 protein subunits, including COPB2.

To place an order, please [Click HERE](#).

Catalog #:A7036

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9276

Isotype:IgG

Swiss Prot:P35606

Purity:Affinity purification

For research use only.