

QKI

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 38kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human QKI

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

QK; Hqk; QK1; QK3; hqkl;

Catalog #: A7043

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9444

Isotype: IgG

Swiss Prot: Q96PU8

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is an RNA-binding protein that regulates pre-mRNA splicing, export of mRNAs from the nucleus, protein translation, and mRNA stability. The encoded protein is involved in myelination and oligodendrocyte differentiation and may play a role in schizophrenia. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).