

TOX

Reactivity:Human

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:58kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human TOX

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TOX1;

Catalog #:A7050

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9760

Isotype:IgG

Swiss Prot:O94900

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene contains a HMG box DNA binding domain. HMG boxes are found in many eukaryotic proteins involved in chromatin assembly, transcription and replication.

This protein may function to regulate T-cell development.

To place an order, please [Click HERE](#).