

GJC1

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:45kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human GJC1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CX45; GJA7;

Catalog #:A7053

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10052

Isotype:IgG

Swiss Prot:P36383

Purity:Affinity purification

For research use only.

Background:

This gene is a member of the connexin gene family. The encoded protein is a component of gap junctions, which are composed of arrays of intercellular channels that provide a route for the diffusion of low molecular weight materials from cell to cell. Alternatively spliced transcript variants encoding the same isoform have been described.

To place an order, please [Click HERE](#).