

NUTF2

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:14kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human NUTF2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NTF2; PP15;

Catalog #:A7057

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10204

Isotype:IgG

Swiss Prot:P61970

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a cytosolic factor that facilitates protein transport into the nucleus. It interacts with the nuclear pore complex glycoprotein p62. This encoded protein acts at a relative late stage of nuclear protein import, subsequent to the initial docking of nuclear import ligand at the nuclear envelope. It is thought to be part of a multicomponent system of cytosolic factors that assemble at the pore complex during nuclear import.

To place an order, please [Click HERE](#).