

## PIAS3

---

**Reactivity:**Human Rat

**Tested applications:**WB IHC IF

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:100 IF 1:50 - 1:100

**Calculated MW:**68kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human PIAS3

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

ZMIZ5;

**Background:**

This gene encodes a member of the PIAS

*To place an order, please [Click HERE](#).*

**Catalog #:**A7060

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**10401

**Isotype:**IgG

**Swiss Prot:**Q9Y6X2

**Purity:**Affinity purification

For research use only.