

VTI1B

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 27kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human VTI1B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

VTI1; VTI2; VTI1L; v-SNARE; vti1-rp1; VTI1-LIKE;

To place an order, please [Click HERE](#).

Catalog #: A7062

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10490

Isotype: IgG

Swiss Prot: Q9UEU0

Purity: Affinity purification

For research use only.