

## AGR2

---

**Reactivity:**Human

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**20kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human AGR2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

AG2; GOB-4; HAG-2; XAG-2; PDIA17; HEL-S-116;

*To place an order, please [Click HERE](#).*

**Catalog #:**A7064

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**10551

**Isotype:**IgG

**Swiss Prot:**O95994

**Purity:**Affinity purification

For research use only.