

FASTK

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:61kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human FASTK

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FAST;

Catalog #:A7069

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10922

Isotype:IgG

Swiss Prot:Q14296

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the serine/threonine protein kinase family. This kinase was shown to be activated rapidly during Fas-mediated apoptosis in Jurkat cells. In response to Fas receptor ligation, it phosphorylates TIA1, an apoptosis-promoting nuclear RNA-binding protein. The encoded protein is a strong inducer of lymphocyte apoptosis. Two transcript variants encoding different isoforms have been found for this gene. Other variants exist, but their full-length natures have not yet been determined.

To place an order, please [Click HERE](#).