

TUSC2

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:12kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human TUSC2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PAP; FUS1; PDAP2; C3orf11;

Background:

This gene is a highly conserved lung cancer candidate gene. No other information about this gene is currently available.

To place an order, please [Click HERE](#).

Catalog #:A7075

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:11334

Isotype:IgG

Swiss Prot:O75896

Purity:Affinity purification

For research use only.