

P2RX2

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 52kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human P2RX2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

P2X2; DFNA41;

Catalog #: A7078

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 22953

Isotype: IgG

Swiss Prot: Q9UBL9

Purity: Affinity purification

For research use only.

Background:

The product of this gene belongs to the family of purinoceptors for ATP. This receptor functions as a ligand-gated ion channel. Binding to ATP mediates synaptic transmission between neurons and from neurons to smooth muscle. Multiple transcript variants encoding distinct isoforms have been identified for this gene.

To place an order, please [Click HERE](#).