

ICOSLG

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:33kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human ICOSLG

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

B7H2; GL50; B7-H2; B7RP1; CD275; ICOSL; LICOS; B7RP-1; ICOS-L;

To place an order, please [Click HERE](#).

Catalog #:A7080

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:23308

Isotype:IgG

Swiss Prot:O75144

Purity:Affinity purification

For research use only.