

FKBP8

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:45kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human FKBP8

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bf

Synonym:

FKBP38; FKBP38;

Catalog #:A7085

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:23770

Isotype:IgG

Swiss Prot:Q14318

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the immunophilin protein family, which play a role in immunoregulation and basic cellular processes involving protein folding and trafficking.

Unlike the other members of the family, this encoded protein does not seem to have

PPlase/rotamase activity. It may have a role in neurons associated with memory function.

To place an order, please [Click HERE](#).