

CLIC4

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:29kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human CLIC4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

H1; huH1; p64H1; CLIC4L; MTCLIC;

Catalog #:A7088

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:25932

Isotype:IgG

Swiss Prot:Q9Y696

Purity:Affinity purification

For research use only.

Background:

Chloride channels are a diverse group of proteins that regulate fundamental cellular processes including stabilization of cell membrane potential, transepithelial transport, maintenance of intracellular pH, and regulation of cell volume. Chloride intracellular channel 4 (CLIC4) protein, encoded by the CLIC4 gene, is a member of the p64 family; the gene is expressed in many tissues and exhibits a intracellular vesicular pattern in Panc-1 cells (pancreatic cancer cells).

To place an order, please [Click HERE](#).