

DNM3

Reactivity:Human Mouse Rat

Tested applications:WB IHC IP

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IP 1:20 - 1:50

Calculated MW:98kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human DNM3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

Dyna III;

Background:

This gene encodes a member of a family of guanosine triphosphate (GTP)-binding proteins that associate with microtubules and are involved in vesicular transport. The encoded protein functions in the development of megakaryocytes. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).

Catalog #:A7092

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:26052

Isotype:IgG

Swiss Prot:Q9UQ16

Purity:Affinity purification

For research use only.