

LDLRAP1

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 34kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human LDLRAP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ARH; ARH1; ARH2; FHCB1; FHCB2;

Catalog #: A7093

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 26119

Isotype: IgG

Swiss Prot: Q5SW96

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a cytosolic protein which contains a phosphotyrosine binding (PTD) domain. The PTD domain has been found to interact with the cytoplasmic tail of the LDL receptor. Mutations in this gene lead to LDL receptor malfunction and cause the disorder autosomal recessive hypercholesterolaemia.

To place an order, please [Click HERE](#).