

TRMT2A

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 69kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human TRMT2A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

fj

Synonym:

HTF9C;

Catalog #: A7095

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 27037

Isotype: IgG

Swiss Prot: Q8IZ69

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is of unknown function. However, it is orthologous to the mouse Trmt2a gene and contains an RNA methyltransferase domain. Expression of this gene varies during the cell cycle, with aberrant expression being a possible biomarker in certain breast cancers. Several transcript variants encoding two different isoforms have been found for this gene.

To place an order, please [Click HERE](#).