

## TRMT2A

---

**Reactivity:**Human

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**69kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic peptide of human TRMT2A

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

fj

**Synonym:**

HTF9C;

**Catalog #:**A7095

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**27037

**Isotype:**IgG

**Swiss Prot:**Q8IZ69

**Purity:**Affinity purification

For research use only.

**Background:**

The protein encoded by this gene is of unknown function. However, it is orthologous to the mouse Trmt2a gene and contains an RNA methyltransferase domain. Expression of this gene varies during the cell cycle, with aberrant expression being a possible biomarker in certain breast cancers. Several transcript variants encoding two different isoforms have been found for this gene.

**To place an order, please [Click HERE](#).**