

SH3GLB1

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:41kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human SH3GLB1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bch

Synonym:

Bif-1; CGI-61; PPP1R70; dJ612B15.2;

Catalog #:A7104

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:51100

Isotype:IgG

Swiss Prot:Q9Y371

Purity:Affinity purification

For research use only.

Background:

This gene encodes a SRC homology 3 domain-containing protein. The encoded protein interacts with the proapoptotic member of the Bcl-2 family, Bcl-2-associated X protein (Bax) and may be involved in regulating apoptotic signaling pathways. This protein may also be involved in maintaining mitochondrial morphology. Alternate splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).