

## SH3GLB1

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC IF

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

**Calculated MW:** 41kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human SH3GLB1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

bch

**Synonym:**

Bif-1; CGI-61; PPP1R70; dJ612B15.2;

**Catalog #:** A7104

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 51100

**Isotype:** IgG

**Swiss Prot:** Q9Y371

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes a SRC homology 3 domain-containing protein. The encoded protein interacts with the proapoptotic member of the Bcl-2 family, Bcl-2-associated X protein (Bax) and may be involved in regulating apoptotic signaling pathways. This protein may also be involved in maintaining mitochondrial morphology. Alternate splicing results in multiple transcript variants.

*To place an order, please [Click HERE](#).*