

SERPINA10

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100 IF 1:50 - 1:100

Calculated MW: 51kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SERPINA10

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PZI; ZPI;

Catalog #: A7106

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 51156

Isotype: IgG

Swiss Prot: Q9UK55

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene belongs to the serpin family. It is predominantly expressed in the liver and secreted in plasma. It inhibits the activity of coagulation factors Xa and XIa in the presence of protein Z, calcium and phospholipid. Mutations in this gene are associated with venous thrombosis. Alternatively spliced transcript variants have been found for this gene.

To place an order, please [Click HERE](#).