

NDE1

Reactivity:Human Mouse

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:39kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human NDE1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1 mg/ml

Synonym:

NDE; LIS4; MHAC; NUDE; NUDE1; HOM-TES-87;

Catalog #:A7112

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:54820

Isotype:IgG

Swiss Prot:Q9NXR1

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the nuclear distribution E (NudE) family of proteins. The encoded protein is localized at the centrosome and interacts with other centrosome components as part of a multiprotein complex that regulates dynein function. This protein plays an essential role in microtubule organization, mitosis and neuronal migration. Mutations in this gene cause lissencephaly 4, a disorder characterized by lissencephaly, severe brain atrophy, microcephaly, and severe mental retardation. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).