

## PTCRA

---

**Reactivity:**Human

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**29kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human PTCRA

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

PTA; PT-ALPHA;

**Catalog #:**A7150

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**171558

**Isotype:**IgG

**Swiss Prot:**Q6ISU1

**Purity:**Affinity purification

For research use only.

**Background:**

The protein encoded by this gene is a single-pass type I membrane protein that is found in immature but not mature T-cells. Along with TCRB and CD3 complex, the encoded protein forms the pre-T-cell receptor complex, which regulates early T-cell development. Four transcript variants encoding different isoforms have been found for this gene.

*To place an order, please [Click HERE](#).*