

CNPY3

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:31kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human CNPY3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CAG4A; ERDA5; TNRC5; PRAT4A;

Background:

PRAT4A is associated with the immature form of TLR4 (MIM 603030) and regulates its cell surface expression .

To place an order, please [Click HERE](#).

Catalog #:A7176

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10695

Isotype:IgG

Swiss Prot:Q9BT09

Purity:Affinity purification

For research use only.