

ELMO3

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 81kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ELMO3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CED12; CED-12; ELMO-3;

Catalog #: A7182

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 79767

Isotype: IgG

Swiss Prot: Q96BJ8

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is similar to a *C. elegans* protein that functions in phagocytosis of apoptotic cells and in cell migration. Other members of this small family of engulfment and cell motility (ELMO) proteins have been shown to interact with the dedicator of cyto-kinesis 1 protein to promote phagocytosis and effect cell shape changes.

To place an order, please [Click HERE](#).