

ETV1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 55kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ETV1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ER81;

Catalog #: A7184

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2115

Isotype: IgG

Swiss Prot: P50549

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the ETS (E twenty-six) family of transcription factors. The ETS proteins regulate many target genes that modulate biological processes like cell growth, angiogenesis, migration, proliferation and differentiation. All ETS proteins contain an ETS DNA-binding domain that binds to DNA sequences containing the consensus 5'-CGGA

To place an order, please [Click HERE](#).